

Modified Approach to Unknowns

Resource Type: Curriculum: Laboratory

Publication Date: 2/20/2003

Authors

Ethel Matthews
Midland College
Midland, TX
Email: ematthews@midland.edu

Paul Mangum
Midland College
Midland, TX

Denise Humphrey
Midland College
Midland, TX

Abstract

This exercise has students apply skills and techniques learned in the lab to identify assigned, but unnamed, bacteria in a defined time frame with emphasis on wise use of media and supplies. This method reduces student involvement in the actual manipulation of microbes and media but allows the opportunity for the student to grasp the process of just how one approaches the identification of a microbe. Traditionally students are assigned unknown microbes to identify with generous media and other supplies to use and liberal access to the microbiology laboratory. The approach described here represents a substantial savings in supplies and considers availability of the microbiology lab. Some colleges have a number of different kinds of lab classes scheduled for one lab room (not only microbiology, but also general biology, anatomy and physiology, genetics, etc.) thus restricting access to the facility outside designated class hours.

Activity

Invitation for User Feedback. If you have used the activity and would like to provide feedback, please send an e-mail to MicrobeLibrary@asmusa.org. Feedback can include ideas which complement the activity and new approaches for implementing the activity. Your comments will be added to the activity under a separate section labeled "Feedback." Comments may be edited.

INTRODUCTION

Background.

This exercise takes place after students have done Gram stains and set up and read primary media and biochemicals. Information regarding the organisms, media use, charts, etc., has been presented in exercises specifically dedicated to the *Enterobacteriaceae* and glucose nonfermenters, *Staphylococci*, *Streptococci*, and related organisms. Students have used various identification charts for all these categories of organisms and have been encouraged to keep accurate notes for all media and supplies.

PROCEDURE

Materials.

- Gram stain reagents, microscope slides, hydrogen peroxide for each student.
- Stock cultures of *Staphylococci*, *Streptococci*, *Enterobacteriaceae*, *Pseudomonas*, *Acinetobacter*.
- Tryptic soy agar (TSA) plates and MacConkey agar plates.
- Blood agar plates, biochemicals, disks, and reagents so that there is one full identification set up for each organism used for unknowns.
- [Name that Bug!](#) Directions and worksheet for students
- [Media Order List \(HTML\)](#)
- [Media Order List \(PDF\)](#)
- [Observation Forms for Day Two \(HTML\)](#)
- [Observation Forms for Day Two \(PDF\)](#)
- [Flow Chart to Identify Streptococci \(HTML\)](#)
- [Flow Chart to Identify Staphylococci \(HTML\)](#)
- [Chart for Identification of Some Enterobacteriaceae and Some Glucose Nonfermenting Gram-negative Bacilli \(HTML\)](#)

Student Version.

Day 1:

Each student receives four pure cultures of bacteria. A *Streptococcus* and a *Staphylococcus* have been streaked on tryptic soy agar (TSA) plates. A gram-negative bacilli from the *Enterobacteriaceae* and/or a nonglucose fermenter such as *Pseudomonas* or *Acinetobacter* have been streaked on MacConkey agar plates.

Determine Gram stain and, if necessary, catalase results. Since TSA is not a selective medium, organisms on these plates must be Gram stained. Recall the use of MacConkey agar to record the Gram stain result without actually performing the Gram stain. You may perform a Gram stain on the MacConkey microbes if you choose. At this point you are expected to submit a list of media appropriate to identifying each organism. The media list is in chart form with areas to designate whether the medium is a plate, broth, slant, or deep; whether it should be streaked confluent, stabbed, fishtailed on the surface, etc.; and whether any additional reagents such as Kovács reagent, disks, etc., are needed. This concludes lab activities for Day 1.

The organisms used in the unknowns are those that have been introduced in previous labs and can be identified using techniques and procedures from previous labs. This is an open book activity and you may use your lab manuals, notes, and identification charts; however, you may not borrow someone else's materials or involve another student in the process.

Day 2:

Visit only the workstations where each of your four organisms have been set up. Fill out a chart that includes reactions (change in media color, hemolysis on blood agar, zone sizes, etc.), interpretation of the results (positive or negative, susceptible or resistant, etc.), and the scientific name of the organism.

Note: samples of the directions and charts given to the students follow in the Supplementary Materials section.

Instructor Version.

This exercise takes place during two 2-hour labs. Students work independently but may use their notes and lab manual.

The instructor assigns a number to each microbe and streaks two microbes to either TSA or MacConkey agar the day before the exercise begins. Note: the objective of Day 1 is to have the students determine Gram stain results and catalase reactions. It is acknowledged that not all of the *Streptococci* will grow on TSA but to streamline set up time and conserve media the instructor may set up the "Streps" using *Enterococcus* species to simulate all the Streps. *Enterococci* will grow on TSA, produce gram-positive cocci on the Gram stain, and show a negative catalase reaction. If the instructor feels that it is more appropriate, the Streps can be set up on blood agar plates, and students can read hemolysis on Day 1. The grading process would need to be adjusted to reflect this change. Schedule at least an hour to do the set up.

Suggested microorganisms:

Streptococcus pyogenes
Escherichia coli
Streptococcus agalactiae
Klebsiella pneumoniae
Streptococcus viridans
Enterobacter aerogenes
Streptococcus pneumoniae
Serratia marcescens
Enterococcus faecalis
Proteus vulgaris
Staphylococcus aureus
Salmonella species
Staphylococcus epidermidis
Shigella flexneri
Staphylococcus saprophyticus
Providentia startii
Pseudomonas aeruginosa
Acinetobacter anitratus

You may want to use a couple of these more than once.

Work stations are often set up with four students per table. Here is an easy way to be certain that each student at the table receives a unique set of organisms. Number the Streps 1 through 6, the Staphs 7 through 10, and the gram-negative bacilli species 11 through 20. Assign each student four organisms as follows:

Table 1	Table 2	Table 3	Table 4	Table 5
1, 7, 11, 16	5, 7, 15, 20	3, 7, 14, 20	1, 8, 13, 20	5, 8, 12, 20
2, 8, 12, 17	6, 8, 11, 17	4, 8, 15, 16	2, 9, 14, 16	6, 9, 13, 16
3, 9, 13, 18	7, 9, 12, 18	5, 9, 11, 18	3, 10, 15, 17	1, 10, 14, 18
4, 10, 14, 19	2, 10, 13, 19	6, 10, 12, 19	4, 7, 11, 19	2, 7, 15, 19

The instructor sets up appropriate media to identify each organism of the entire set of organisms that were provided to the students (about 20 in all) prior to Day 2. It will take between 2 and 3 hours to label and inoculate the plates and tubes, place disks, etc., provided that the media is available. Before students arrive on Day 2, all media for organism #1 will be set at workstation #1, all media for organism #2 will be set at workstation #2, etc., until all media is at the individual stations. The instructor adds Kovács and ferric chloride reagent to the sulfur-indole motility (SIM) and phenylalanine tubes, respectively, as needed, just prior to the arrival of the students. It is advisable to set up a SIM and phenylalanine for each class period for each of the gram negatives as the colors of the added reagents fade or absorb into the media after awhile.

Safety Issues.

Students should be reminded that they are working with bacteria that are either recognized pathogens or have the potential to be pathogens. Students should observe aseptic technique and safe handling procedures for the microbes, safe use of the

Bunsen burner, and concern for flammable reagents such as the acetone-alcohol decolorizer for Gram stain decolorization.

ASSESSMENT and OUTCOMES

Suggestions for Assessment.

The student can receive a maximum of 25 points for identifying each organism, based on Gram stain and catalase results (5 points), appropriate order for media (10 points), accurate observation and interpretation of media that is set up (10 points), and accurate identification of microbe (5 points). The grade for ordering media is broken down into: media required, proper set-up method, and additional reagents and disks needed. The grade for observation and interpretation is broken down to each biochemical, etc. This method allows for an objective and relatively easy grading approach to unknowns.

If students accurately order appropriate media (including proper notations) for Day 1, they receive all points. An inadequate request results in loss of the points for the missing media. Ordering too much media is penalized by half points.

Field Testing.

This approach to unknowns has been used for both regular and summer sessions for the past 7 years. Class sizes vary from 16 to 24 students. Student comments at the end of the exercise indicate that they finally see the big picture, and a few have said that it was their favorite exercise of the semester.

Student Data.

Example of a student comment:

"Unknowns - Testing my ability to see if I really understood how to identify a bug - putting it altogether....It was good to see that I could critically think a bit and come up with the right bug!"

[Sample of Student Outcomes](#)

Individual GIF images of samples of student outcomes:

- [Student sample of Name That Bug worksheet](#)
- [Student sample 2 of Media Order List worksheet](#)
- [Student sample 3 of Media Order List worksheet](#)
- [Student sample 4 of Media Order List worksheet](#)
- [Student sample of Gram Positive worksheet](#)
- [Student sample of Gram Negative worksheet](#)

SUPPLEMENTARY MATERIALS

Answer Keys.

- [Key for Gram-negative Bacilli \(HTML\)](#)
- [Key for Gram-positive Cocci - Catalase Negative \(HTML\)](#)
- [Key for Gram-positive Cocci - Catalase Positive\(HTML\)](#)

Curriculum Resources

Modified Approach to Unknowns - Matthews

Name That Bug!

Organism #	Gram stain report	Catalase (if applicable)

Directions:

This is a test of your ability to identify microorganisms. You may use your notes, handouts, lab manual, etc., but you must work independently. You **may not** consult with your neighbor, use his or her notes, ask his opinion, or ask someone to focus your microscope. Assume that I am the patient and I will become very agitated if you expect me to know what is wrong with me. YOU are the professional, you figure it out.

Day 1:

Today you need to determine the Gram stain of your organism. Remember, let the media do some of the work for you. On **gram positive** organisms, you need to perform a catalase test. If you do a catalase test on **gram negative** organisms, please be prepared to explain how that information is relevant to the identification of the organism. Using the chart provided, please **circle** the agar and/or test you are requesting. **Indicate** whether this agar is a broth, slant, or plate. **Indicate** the tool you will be using to inoculate this medium, and the method of inoculation (isolation, confluent, stab, etc). **Check** all other things that apply such as pre- or post-incubation additives are required. **Describe** test set up on the appropriate line or at the page bottom for all asterisked (*) tests, such as CAMP, coagulase, etc. **Assume that you have 24 hours** to identify your organisms. Your list should be sufficient for identification without undue waste of materials and time. Turn in ALL worksheets with Gram stain results, catalase testing, and the media list **at the end of this class period.**

Day 2:

Next class period. The appropriate work will have been set up in order for you to identify your organisms. Regardless of what you may have requested, you can still identify your organisms to species level (with few exceptions). You can proceed to the workstation that corresponds to your organism number and read and interpret **all** tests or media you find present. Record on your worksheets both the **observations** (what you actually saw) and the interpretations (what you think this means: positive or negative (+/-), susceptible or resistant (S/R), etc.). **Include** all disk measurements. In all final identifications, the organism names should be correctly spelled with genus and species underlined by you.

Scoring:

- 5 pts. Gram stain results and catalase test results
- 5 pts. list of media and tests, etc., to identify organism

- 10 pts. reading and interpretation of media, etc.
- 5 pts. organism name

(25 points per organism x 4 organisms = 100 points)

Name _____

Date ____ Class Day/ Time

Organism # _____

Gram Stain

Medium	broth, plate, slant, deep, other	needle, loop, other	Method of inoculation						Preincubation reagents/ disks				Post incubation reagents	
--------	----------------------------------	---------------------	-----------------------	--	--	--	--	--	-------------------------------	--	--	--	--------------------------	--

			Stab	fishtail	isolation	confluent	simple	AP@	AA@	Nov.	<i>Staph aureus</i>	Mineral Oil	Kovac	Fer chl
BAP														
BEA														
CAMP*														
Coag.*														
DNA														
Lysine														
MSA														
TSA														
Phen Deam														
SIM														
Sim.Cit.														
TSI														
Urea														
MacC														
Oxidase*														

* describe in detail

Name _____ Date _____ Class Day /Time _____

Gram Positive Organisms

Report all reactions at the work station using correct and appropriate format (hemolysis, zone size, color, +/-, S/R)

Organism

Blood agar plate hemolysis _____ (Alpha, beta, gamma)

Bacitracin disk ____ (S/ R) Optochin disk ____ (S/ R) Bile esculin agar ____ (+/ -)
Zone size ____ mm Zone size ____ mm Color

CAMP test ____ (+/ -)

Coagulase ____ (+/ -) Mannitol Salt Agar ____ (+/ -) Novobiocin disk ____ (S/ R)
Describe _____ _____ Zone size ____ mm
(Growth on plate/ color of agar)

Organism identification

Organism

Blood agar plate hemolysis _____ (Alpha, beta, gamma)

Bacitracin disk ____ (S/ R) Optochin disk ____ (S/ R) Bile esculin agar ____ (+/ -)
Zone size ____ mm Zone size ____ mm Color

CAMP test ____ (+/ -)

Coagulase ____ (+/ -) Mannitol Salt Agar ____ (+/ -) Novobiocin disk ____ (S/ R)
Describe _____ _____ Zone size ____ mm
(Growth on plate/ color of agar)

Organism identification _____

Gram Negative Organisms

Report all reactions at the work station using correct and appropriate format. Refer to your notes for TSI reactions. Use the first line under the biochemical name to record the interpretation of the reactions and the second line to record your observation (color change, turbidity, etc.) on which you based your interpretation.

Organism

	TSI	SIM	Simmons Citrate	Lysine Decarb.	Urease	Phenyl. Deam	DNase	Oxidase
Interp.	_____	_____	_____	_____	_____			
Observ.	_____	S_____	_____	_____	_____			
		I						
		M						

Organism Identification

Organism


	TSI	SIM	Simmons Citrate	Lysine Decarb.	Urease	Phenyl. Deam	DNase	Oxidase
Interp.	_____	_____	_____	_____	_____			
Observ.	_____	S_____	_____	_____	_____			
		I						
		M						

Organism Identification

Curriculum Resources

Modified Approach to Unknowns - Matthews


Flow Chart to Identify Streptococci



Curriculum Resources

Modified Approach to Unknowns - Matthews

Flow Chart for the Identification of Staphylococci



 Curriculum Resources

Modified Approach to Unknowns - Matthews

**Chart for Identification of Some *Enterobacteriaceae* and
Some Glucose Non-fermenting Gram-negative Bacteria**

organism	TSI	SIM	Citrate	Lysine	Urea	PDA	DNA	Oxidase
<i>Citrobacter freundii</i>	K/A H ₂ S+	+++	+	-	-	-	-	-
<i>Enterobacter aerogenes</i>	A/A +	+++	+	+	-	-	-	-
<i>Escherichia coli</i>	A/A +	+++	-	±	-	-	-	-
<i>Klebsiella pneumoniae</i>	A/A	+++	+	+	±	-	-	-
<i>Morganella morganii</i>	K/A H ₂ S+	+++	-	-	+	+	-	-
<i>Proteus vulgaris</i>	K/A H ₂ S+	+++	-	-	+	+	-	-
<i>Providentia stuartii</i>	K/A	++±	+	-	-	+	-	-
<i>Salmonella</i> species	K/A H ₂ S+	+++	+	+	-	-	-	-
<i>Serratia marcescens</i>	K/A	++	+	+	-	-	+	-
<i>Shigella flexneri</i>	K/A	---	-	-	-	-	-	-
<i>Pseudomonas aeruginosa</i> *	K/K	++	+	-	±	-	-	+
<i>Acinetobacter anitratus</i> **	K/K	---	-	-	±	-	-	-

NOTES

None of the *Enterobacteriaceae* are oxidase positive.

Occasionally some of the *Enterobacteriaceae* that usually produce A/A results on the TSI slants may show a slight amount of alkaline reaction on the slant. Read as A/A.

Glucose nonfermenters:

* *Pseudomonas aeruginosa* is beta-hemolytic on blood agar plates, oxidase positive, and often produces a green pigment. Motility is difficult to read as this organism is aerophilic.

** *Acinetobacter anitratus* is nonhemolytic on blood agar plates, oxidase negative, and does not produce

pigment.

Abbreviations: TSI=triple sugar iron, SIM=sulfur-indole motility, Citrate=Simmons citrate, Lysine=lysine decarboxylase, Urea=Christensen=s urea, PDA= phenylalanine deaminase, DNA= DNase.

Note to Instructors

Biochemical reactions of some of these organisms could actually be designated as v or + according to references such as *Bailey and Scott's Diagnostic Microbiology*, however, the strains in use here at Midland College are giving the reactions as designated in the above chart. The goal of the unknowns exercise is to get the student to use an identification chart with a limited number of biochemical reactions. Also having a number of organisms included in the chart, even though not all are used in the assigned unknowns, is a way of expanding the critical thinking skills of the student.

Curriculum Resources

Modified Approach to Unknowns - Matthews Sample of Student Outcomes

Name



97

NAME THAT BUG!

Organism #	Gram Stain Reaction	Catalase (if applicable)
20	gram positive ^{negative} bacilli	—
17.2 ^{list}	gram positive ^{negative} bacilli	—
37.2	gram positive cocci	—
60	Gram positive cocci	+
97.2		

Directions:

This is a test of your ability to identify microorganisms. You may use your notes, handouts, lab book, etc... but you must work independently. You **MAY NOT** consult with your neighbor, use his/her notes, ask their opinion or ask someone to focus your microscope. Assume that I am the patient and I will become very agitated if you expect me to know what is wrong with me. **YOU** are the professional, you figure it out.

DAY 1: Today. You will need to determine the gram stain of your organism.

Remember, let the media do some of the work for you. On **gram positive** organisms you will need to perform a catalase test. If you do a catalase test on **gram negative** organisms please be prepared to **explain how that information is relevant to the identification of the organism**. Using the chart provided, please circle the agar/test you are requesting. **Write in** whether this agar is a broth, slant, or plate. **Write in** what tool you will be using to inoculate this media. **Check** the method of inoculation, isolate, confluent, etc... **Check ALL** other things that apply, if pre or post incubation additives are required. **Describe test set-up** in the appropriate line or at the bottom for all asterick (*) tests, such as CAMP, Coagulase, etc... **Assume you have 24 hours** to identify your organisms. Your list should be sufficient without undue waste of materials and time. Turn in ALL worksheets with gram stain results, catalase testing, and media list **at the end of this class period**.

DAY 2: Next class period. The appropriate work will have been set up in order for you to identify your organisms. Regardless of what you may have requested, you can still identify your organisms to species level, with rare exceptions. You can proceed to the work station that corresponds to your organism number and read/interpret **all** tests or media you find present. **Record** on your worksheets both the **observations** (what you actually saw) and the **interpretation** (what you think this means, +/-, S/R, etc...). **Include** all disk measurements. All final identifications/ organism names should be correctly spelled with genus and species underlined by you.

Scoring: 5 pts -- gram stain reaction / catalase
 5 pts -- list of media / tests, etc... to identify organism
 10 pts -- reading & interpretation of media/ biochemical, etc..
 5 pts -- organism name
 (25 points per organism X 4 = 100 points) _____

NAME [Redacted] DATE 10/31/01 CLASS DAY/TIME MW 8:00
 ORGANISM # 2 GRAM STAIN Gram positive cocci

Agar	Method of Inoculation			Pre-incubation			Post-incubation REAGENTS						
	Broth, plate Or slant	Needle Or Loop	Stab	Isolation	confluent	Simple/Fairhill	"P"	"A"	Nov.	Staph. aureus	Mineral oil	Ferric Chloride	Kovac
BAP	Plate	loop			confluent		P	A					
BEA	Plate	loop			confluent								
CAMP	Plate	loop			Streak straight					Staph aureus			
COAGULASE					confluent								
DINA													
LYSINE													
MSA													
TSA													
Phen. Dec.													
SIM													
Simon's Citrate													
TSI													
UREA													
MAC													
OXIDASE*													

+C
+4

10/31/01

* Describe set-up in detail.
 on BAP Streak Staph aureus
 Streak group A streaked at 90°
 $\frac{10}{14} = 0.714 = 3.6$

NAME  DATE 12/31/01 CLASS DAY/TIME MW 8:00

ORGANISM # 8 GRAM STAIN gram positive cocci

Agar	Method of Inoculation			Pre-Incubation			Post-Incubation REAGENTS	
	Broth, plate Or slant	Needle Or Loop	Streak	DISKS	Staph. aureus	Mineral oil	Ferric Chloride	Kovac
BAP	plate	loop	confluent	*P*	Nov.			
BEA			Non-fluor		Nov.			
CAMP*								
COAGULASE	plate	loop						
DNA								
LYSINE								
MISA	plate	loop						
TSA								
Phen. Dea.								
SIM								
Simon's Citrate								
TSI								
UREA								
MAC								
OXIDASE*								

* Describe set-up in detail.

Handwritten notes and signatures

NAME: [REDACTED] DATE: 10/31/01 CLASS DAY/TIME: MW 8:00
 ORGANISM # 12 GRAM STAIN: gram negative bacilli

Agar	Method of Inoculation				Pre-incubation			Post-incubation					
	Broth, plate Or slant	Needle Or Loop	Stab	Isolation confluent	Streak	Simple/Fish tail	"P"	"A"	Nov.	Staph. aureus	Mineral oil	Ferric Chloride	Kovac
BAP													
BEA													
CAMP													
COAGULASE													
DNA	plate	needle				Simple							
LYSINE	broth	needle											
MSA													
TSA													
Phen. Dea.													
SIM	slant	needle	✓			Fish tail							
Simon's Citrate	deep	needle	✓										
TSI	slant	needle	✓			Simple							
UREA	slant	needle	✓			Fish tail							
MAC	slant	needle				Fish tail							
OXIDASE*	plate	loop				Fish tail							
	ON	com											
	stab												

This substance #12 + #20 per gram negative bacilli - both lists were identical

Describe set-up in detail:
 plate King on Swab
 drop oxidase reagent
 Purple(+) positive
 NCC (-) negative
 $\frac{35}{37} \approx 4.7$

Date 11/5 Class day/time 8:00AM**GRAM POSITIVE ORGANISMS**

Report all reactions at the work station using correct and appropriate Format. (zone size, color, positive (+), negative (-), susceptible, resistant, etc...

ORGANISM # 2Blood plate agar hemolysis beta (alpha, beta, or gamma)

Bacitracin disk S (S/R) Optochin disk R (S/R) Bile Esculin Agar - (+/-)
 Zone size 13 mm Zone size 2 mm ~~zone~~ size 100 mm

CAMP Test - (+/-)color
of mediaCoagulase - (+/-)
(Acid-staph)Mannitol Salt Agar -Novobiocin disk - (S/R)
Zone size - mm

(growth on plate? Color of agar?)

Organism identification group A Beta Streptococcus

15

ORGANISM # BBlood plate agar hemolysis beta (alpha, beta, or gamma)

Bacitracin disk - (S/R) Optochin disk - (S/R) Bile Esculin Agar - (+/-)
 Zone size - mm Zone size - mm color -

CAMP Test - (+/-)Coagulase + (+/-)Mannitol Salt Agar +Novobiocin disk S (S/R)
Zone size 18 mmyes, yellow
(growth on plate? Color of agar?)Organism identification Staphylococcus aureus

15

GRAM NEGATIVE ORGANISMS

Report all reactions at the work station using correct and appropriate format, such as for TSI reactions (refer to your notes). Use the first line under the biochemical name to record the interpretation of your reactions and the second line to record your observation (color change, turbidity, etc...) on which you based your interpretation.

ORGANISM # 20

	TSI	SIM	Simmons Citrate	Lysine Decarb	Urease	Phenyl. Deam.	Dnase	Oxidase
Intrep. K/A H ₂ S	++	+	+	+	-	=	=	-
Observ.	S black	blue	blue	purple	NCC	NCC	no zone of clearing	-
	I NCC							
	M turbid							

15 ORGANISM IDENTIFICATION Salmonella species

ORGANISM # 12

	TSI	SIM	Simmons Citrate	Lysine Decarb	Urease	Phenyl. Deam.	Dnase	Oxidase
Intrep. K/A	++	+	+	+	-	+	=	-
Observ.	S NCC	blue	blue	yellow	NCC	green	no zone of clearing	-
	I red							
	M turbid							

15 ORGANISM IDENTIFICATION Providencia stuartii

Curriculum Resources

**Modified Approach to Unknowns - Matthews
Media Order Key for Gram-negative Bacilli**

NAME Key DATE _____ CLASS DAY/TIME _____

ORGANISM # _____ GRAM STAIN gram negative bacilli

Agar	Broth, plate Or slant	Needle Or Loop	Method of Inoculation				Pre-incubation				Post-incubation REAGENTS		
			Stab	Isolation	confluent	Simple / Fishtail	*P*	*A*	Nov.	Staph. aureus	Mineral oil	Ferric Chloride	Kovac
BAP													
BEA													
CAMP *													
COAGULASE Accu-staph*													
4 DNA	plate	needle				simple						5-6 drops	
5 LYSINE	broth	needle				simple						5-6 drops	
MSA													
TSA													
5 Phen. Dea.	plate	needle				stab							
5 SIM	deep	needle stab											8-10 drops
4 Simon's Citrate	plate	needle				simple							
5 TSI	plate	needle stab				stab fish tail							
4 UREA	plate	needle				stab fish tail							
MAC													
5 OXIDASE*	broth swab												

* Describe set-up in detail.
 1. touch swab on colony
 2. drop 1 drop of oxidase reagent on swab

Curriculum Resources

Modified Approach to Unknowns - Matthews Media Order Key for Gram-positive Cocci - Catalase Negative

NAME Key DATE _____ CLASS DAY/TIME _____

ORGANISM # _____ GRAM STAIN gram positive cocci (catalase negative)

Agar	Broth, plate Or slant	Needle Or Loop	Method of Inoculation				Pre-incubation				Post-incubation REAGENTS		
			Stab	Isolation	confluent	Simple / Fishtail	*P*	*A*	Nov.	Staph. aureus	Mineral oil	Ferric Chloride	Kovac
<u>BAP</u>	<u>plate</u>	<u>loop</u>			✓		✓	✓					
<u>BEA</u>	<u>plate</u>	<u>needle</u>				<u>simple</u>							
<u>CAMP</u>		<u>loop</u>									✓		
COAGULASE Accu-staph*													
DNA													
LYSINE													
MSA													
TSA													
Phen. Dea.													
SIM													
Simon's Citrate													
TSI													
UREA													
MAC													
OXIDASE*													

BAP

 Unknown
Staph aureus

* Describe set-up in detail.

Curriculum Resources

Modified Approach to Unknowns - Matthews Media Order Key for Gram-positive Cocci - Catalase Positive

NAME Key DATE _____ CLASS DAY/TIME _____

ORGANISM # _____ GRAM STAIN gram positive cocci (catalase positive)

Agar	Method of Inoculation		Pre-Incubation						Post-Incubation REAGENTS				
	Broth, plate Or slant	Needle Or Loop	Stab	Isolation	STREAK	Simple / Fishtail	*P*	*A*	Nov.	Staph. aureus	Mineral oil	Ferric Chloride	Kovac
5 BAP	plate	loop			✓				✓				
BEA													
CAMP *													
5 COAGULASE Accu-staph* BNA	inoculated	loop			place 1 drop of Accu-Staph reagent on well make organism in drop first & loop rock / held up to 1 min.								
LYSINE													
4 MSA	plate	loop			✓								
TSA													
Phen. Dea.													
SIM													
Simon's Citrate													
TSI													
UREA													
MAC													
OXIDASE*													

* Describe set-up in detail.